

AQUASAFE® WSL50 & 25 PLUS

Intermediate portable drinking water laboratories

The Aquasafe® WSL50 Plus is an advanced portable water testing laboratory offering physicochemical and microbiological testing.

The WSL50 Plus is an easy to use portable kit, giving a range of accurate results in a single, easy to transport, hard wearing, waterproof and buoyant case. Contains a powerful digital dual chamber microbiological incubator that docks simply into the case to allow use in both the field or the laboratory. A rugged photometer with testing capabilities for over 40 parameters is also included in the WSL50 Plus. All consumables and accessories are protected by a specifically engineered rack system mounted in the case.

The Aquasafe® WSL25 Plus is a more compact, simple kit which stil offers an excellent range of microbiological and physicochemical tests using familiar field testing methods. Contains a powerful single chamber microbiological incubator, flexible comparator block with wheels and easy to use pocket pH and EC meters.







AquaSafe

WSL50 Plus

W.H.O Guidelines



BENEFITS

- Accurate digital temperature control to 0.1°C ensures optimal bacterial growth
- Process samples while transporting and carrying the portable laboratory with the upright incubator format. Samples remain in the correct orientation even when being carried
- The WSL25 Plus's small lightweight case gives an easy portable solution, even for just one user



- Up to 24 hour countdown timer that pauses in the event of power failure, ensuring no false negative microbiological results
- Total control over incubation temperature.
 Other cultures are possible using different selective media
- PetriLok® petridish carriage system stops incubator contamination from the user's hands and ensures stability of dishes in transit, avoiding media spillage



 Powerful single or double incubator docking station gives flexibility to switch between in-case or on bench use



- Save time with fast and easy calibration in the field
- Single waterproof (IP68) cases makes it easy for one user to carry and floats to avoid loss if dropped in water
- Instant physico-chemical results with the easy to use comparator with colour wheels plus pH and EC pocket meters saves time and money in the field



- Never run out of reagents again. In-built removable accessories tray makes identifying consumable and reagent stock levels simple and avoids the possibility of damage or breakage during transport
- The high performance battery pack built into the case is hidden from potential theft behind a locked hinged panel while being easily accessible using the simple tool provided
- Save time and money in the laboratory. Ideal as either the primary in-field testing method OR as a pre-screening tool to reduce the number of samples centrally tested

IN YOUR KIT

| IN YOUR KIT | WSL50 PLUS | OTY | WSL25 PLUS | OTY | | |
|--|---|-----|------------|-----|--|--|
| Single silver waterproof wheeled case with multiple handle options and incubator | 1 | 1 | | | | |
| docking station | | 1 | | | | |
| Single black lightweight IP68 waterproof case with top handle, incubator docking | | | ✓ | 1 | | |
| station and shoulder strap AquaSafe double incubator | / | 1 | | | | |
| - | V | 1 | | | | |
| AquaSafe Single incubator | | | / | 1 | | |
| PetriLock cassette with cap | √ | 2 | / | 1 | | |
| Aluminium petri dishes | √ | 50 | / | 25 | | |
| Membrane Filtration Unit | / | 1 | / | 1 | | |
| Sterilised water for MSLB preparation | ✓ | 6 | √ √ | 6 | | |
| Sample collection bottle with dechlorination tablets | tion bottle with dechlorination tablets | | | | | |
| Sterile membrane filters (pk100) | ✓ | 3 | / | 3 | | |
| Sterile absorbent pads (pk100) | ✓ | 3 | ✓ | 3 | | |
| Membrane Lauryl sulfate broth sachet | ✓ | 6 | 1 | 6 | | |
| HydroLite HL101 pH meter | 1 | 1 | 1 | 1 | | |
| HydroLite HL102 EC meter | / | 1 | / | 1 | | |
| pH 4.01 buffer solution | / | 1 | / | 1 | | |
| pH 7.00 buffer solution | / | 1 | / | 1 | | |
| EC 1413 uS calibration solution | · · | 1 | / | 1 | | |
| Nitrate reagents (pk100) | 1 | 3 | / | 3 | | |
| Nitrate test tube | 1 | 1 | / | 1 | | |
| Test tube brush | 1 | 1 | / | 1 | | |
| Chlorine free/total reagents (pk100) | 1 | 3 | / | 3 | | |
| Ammonia reagets (pk100) | 1 | 3 | / | 3 | | |
| Nitrite reagents (pk100) | / | 3 | / | 3 | | |
| Digital thermometer | √ | | √ | | | |
| | ✓ | 1 | ✓ ✓ | 1 | | |
| Comparator | | 1 | | 1 | | |
| Comparator cells | / | 4 | / | 4 | | |
| Comparator Chlorine disc | ✓ | 1 | ✓ | 1 | | |
| Comparator Nitrite disc | ✓ | 1 | 1 | 1 | | |
| Comparator Ammonia disc | 1 | 1 | 1 | 1 | | |
| Comparator Nitrate disc | ✓ | 1 | / | 1 | | |
| Arsenic Test Kit (300 tests) | ✓ | 1 | ✓ | 1 | | |
| Turbidity Tube | ✓ | 1 | ✓ | 1 | | |
| Multi-purpose screwdriver | ✓ | 1 | 1 | 1 | | |
| Pad dispenser | ✓ | 1 | ✓ | 1 | | |
| Eye Glass | ✓ | 1 | ✓ | 1 | | |
| Forceps | ✓ | 1 | ✓ | 1 | | |
| 30ml dropper bottle for methanol | 1 | 1 | 1 | 1 | | |
| 5m sterile syringe | 1 | 1 | 1 | 1 | | |
| 5m sampling cable with carabiner | 1 | 1 | 1 | 1 | | |
| Calibration block | 1 | 1 | / | 1 | | |
| Silicone grease | 1 | 1 | / | 1 | | |
| Accessory Boxes | 1 | 3 | | | | |
| Sterilisable work surface | 1 | 1 | | | | |
| Vacuum pump | 1 | 1 | / | 1 | | |
| Vaccum tubing | / | 1 | / | 1 | | |
| 12V 4A lead-acid charger/power supply | / | 1 | / | 1 | | |
| Vehicle charging cable | / | 1 | / | 1 | | |
| IEC mains cable - UK type G | 1 | 1 | / | 1 | | |
| UK to EU adaptor | 1 | 1 | / | 1 | | |
| | | | | - | | |

| DESCRIPTION | PACK SIZE | PART NUMBER |
|----------------------|-----------|-------------|
| Aquasafe® WSL50 Plus | Unit | WSL50PLUS |
| Aquasafe® WSL25 Plus | Unit | WSL25PLUS |



TECHNICAL SPECIFICATIONS

| | ws | L50 PLUS | ws | L25 PLUS | |
|---|----|--|----|----------------------------|--|
| Dimensions (mm) | 0 | 500 x 305 x 457 | 0 | 430 x 290 x 280 | |
| Weight | 0 | 20kg | 0 | 12kg | |
| Approvals | 0 | CE Conformity, IP68 | 1- | | |
| Country of origin | 0 | United Kingdom | | | |
| Digital incubator test programs | 0 | 37°C | | | |
| | 0 | 44°C | | | |
| | 0 | User defined (20°C to 50°C) | | | |
| | 0 | Calibration | | | |
| Timer | 0 | Adjustable 24 hour countdown timer | | | |
| No. incubation chambers | 0 | 2 | 0 | 1 | |
| Maximum petridish capacity | 0 | 50 (25 per incubator chamber) | 0 | 25 | |
| Incubator temperature range | 0 | 20°C to 50°C | | | |
| Incubator temperature resolution | 0 | 0.1°C | | | |
| Incubator temperature stability | 0 | +/- 0.1°C | | | |
| Recommended storage temperature | 0 | 0°C to 70°C | | | |
| Optimum operating temperature | 0 | 0°C to 40°C | | | |
| Under interface | 0 | LCD Display | | | |
| | 0 | Tactile soft keys | | | |
| | 0 | Audible alerts | | | |
| Power supply | 0 | 12-24 V DC | | | |
| | 0 | Removeable lead acid power pack | | | |
| | 0 | (Rechargeable) | | | |
| Power consumption | 0 | 37°C Cycle at 20°C ambient | 0 | 37°C Cycle at 20°C ambient | |
| | 0 | Max power 16.1W | 0 | Max power 8.1W | |
| | 0 | Average power 6W | 0 | Average power 3W | |
| Test parameters included (with reagent) | 0 | Ammonia | | | |
| | 0 | Arsenic (total) | | | |
| | 0 | Chlorine (residual) | | | |
| | 0 | Chlorine (Combined) | | | |
| | 0 | Coliform (TTC or TC) | | | |
| | 0 | Conductivity | | | |
| | 0 | Nitrate | | | |
| | 0 | Nitrite | | | |
| | 0 | рН | | | |
| | 0 | Temperature | | | |
| | 0 | Turbidity | | | |
| Total test parameters possible (with | 0 | Over 50 | | | |
| additional reagent or accessories) | | View to the control of the control o | | 1. | |
| Photometric testing | 0 | Via sturdy comparator block with colour wheels | | | |
| Manufacturer's warranty | 0 | Standard 1 year from purchase | | | |

The Technology Centre, Station Road, Thatcham, Berkshire, United Kingdom, RG19 4HZ